


## **ISWP Static Stability Test**

#### **Assembly Instructions**









# Please read the following document in its entirety before purchasing materials and assembling.

**Design of an ISWP Standards Static Stability** © 2019; University of Pittsburgh.

The International Society of Wheelchair Professionals (ISWP) Static Stability mechanical assembly instructions are made available to the public subject to the following <a href="Creative Commons - Attribution - ShareAlike 4.0 International">Creative Commons - Attribution - ShareAlike 4.0 International</a>. Accordingly, the manual and materials may be downloaded, duplicated, transmitted and otherwise distributed for educational or research purposes, as well as commercially, provided proper credits are given to the University of Pittsburgh and the International Society of Wheelchair Professionals research team. In addition, you must provide a link to the license and also indicate if any changes were made to the materials. If you remix, transform, or build upon the materials, you must distribute your contributions under the same license as the original.

University of Pittsburgh scientists are working with the U.S. Agency for International Development (USAID) under a multi-year sub-award to develop the International Society of Wheelchair Professionals, a global network to ensure a level of standardization, certification and oversight, to teach and professionalize wheelchair services, and to build affiliations to put better equipment in the right hands. Since 2002, USAID has granted more than \$45 million to improve wheelchairs and wheelchair services worldwide. This sub-award – Agreement No. APC-GM-0068 – was presented by Advancing Partners & Communities, a cooperative agreement funded through USAID under Agreement No. AIDOAA-A-12-00047, beginning Oct. 1, 2012.

For further information on use of the ISWP Static Stability assembly instructions, contact the University of Pittsburgh's Innovation Institute at 412-383-7670 or the International Society of Wheelchair Professionals at intlsocietywheelchairprof@gmail.com.









#### **Table of Contents**

Notes:	4
Cutting Notes:	4
Assembly Notes:	4
Tools Required	5
List of drawings	5
Frame	5
Sheet	5
Ramp	5
Stop bar	5
Building the Frame	6
Building the Top Sheet	11
Building the Ramp	12
Building the Stop Bar	13
Bill of Materials	14
Drawings	18









#### **Notes:**

These instructions are to be paired with the dimensioned and assembly drawings for part names and details about each part or assembly.

Assembly instructions may only dictate about one part, but drilling and placement directions apply to all parts of the same name.

All hardware used in this assembly are in ANSI Inch, however ANSI Metric are acceptable alternatives.

All screws used to put lumber together are #8 x3in Screws unless otherwise specified in the Bill of materials.

All dimensions within this set of instructions is given in inches and millimeters.

A part guide has been included at the end of this document.

File down all exposed edges and corners to avoid any nicks or scratches.

#### **Cutting Notes:**

All parts should be cut, and holes drilled before the assembly is started.

#### **Assembly Notes:**

The full assembly is broken down into four subassemblies to simplify explanation, as shown in Figure 1.

Make sure to have all the materials and hardware before assembling.




Figure 1. Four subassemblies. 1. Frame, 2. Sheet, 3. Ramp, 4. Stop bar

For further detail, please see Frame-Bill of Materials and Frame Hardware-Bill of Materials.









#### **Tools Required**

Hammer Philips Screwdriver

Drill

Allen (M8 and M10)

Lumber cutting saw machine

#### **List of drawings**

#### Frame

Part A-Outer Left

Part A-Outer Right

Part A-Inner Left and Inner Right

Part B-Back

Part B-Front

Part C-With rib

Part C

Part D

Part E

Part F

Part G-With mount wheel

Part G-Without mount wheel

#### **Sheet**

Part H

Part I

#### Ramp

Part J

Part K

#### Stop bar

Part L

Part M









#### **Building the Frame**




Figure 2. Frame

See Frame-Bill of Materials and Frame Hardware-Bill of Materials, review each drawing for dimensions before proceeding.

Pay attention to the orientation of the p arts in Figure 3 when building the external part of the frame.




Figure 3. External parts. Frame











Figure 4.Left Side of the Frame

Start by putting together *Part A-Outer Left* and *Part G-Without wheel mount* by screwing them together (Part number 3GCS10BK) on the opposite end of the hole.

Screw Part G and three Parts F to the 3 in casters.

Proceed with screwing *Part G* and all *Parts F* to *Part A-outer Left*, see Figure 4.






Figure 5. Part A- Outer Right









Two *Parts G-without wheel mount* have to be attached to *Part A-Outer Right* on the opposite end of the hole, see Figure 5.

To add the pull handles, use 10-24 x 2.25in Hex Head Screw, 10-24 Locknut, #10 Washers and add them to *Part A-Outer Right* (Figure 5).




Figure 6. Inner ends. Frame

Put Part B-Front and Part B-Back together with Part A-Outer Left and Part A-Outer Right.

Make sure to put Part B-Back at the end where the hole of Part A-Outer Left is.

Add Part A-Inner Left and Right to the frame, see figure 6.










Figure 7. Internal supports. Frame

Add all internal supports (*Parts C*) that go between *Parts A*. Those highlighted in blue are *Parts C-With rib*, see Figure 7.

Finally, add *Parts D* and *Part E*, highlighted in blue in Figure 8. Keep the hole drilled in these parts aligned for the steel rod.




Figure 8. Parts D and Part E. Frame









Add the steel rod and place each shaft collar in the inner edge of the frame (Figure 9).




Figure 9. Steel rod and shaft collars









#### **Building the Top Sheet**




Figure 10. Top sheet

See **BOM Sheet Assembly** before proceeding.

Place *Part H* over the frame and on the opposite end where the steel rod is, using screws (Part number 3GCS10BK).

Similarly, place Part I next to Part H.

Add 3 tiedowns on the long sides of the top sheet and secure them.











#### **Building the Ramp**



Figure 11. Ramp

See **BOM Ramp Assembly** before proceeding.

Make sure to cut *Parts J* and *K* before putting them together using screws (Part number 3GCS10BK).

Finally, add three hinges to hold the ramp and the top sheet together, highlited in Figure 12.




Figure 12. Hinges. Ramp









#### **Building the Stop Bar**




Figure 13. Stop bar

See <u>BOM Stop bar Assembly</u> before proceeding, note this time you need 5 in screws (Part number: 91253A468).

Align *Parts M* and add *Parts L* on each end. Put them together using screws (Part number: 91253A468).

Place the stop bar over the top sheet, see Figure 14.




Figure 14. Stop bar in place









# **Bill of Materials**















### International Society of Wheelchair Professionals I wheelchairnet.ORG 6425 Penn Avenue. Suite 401. Pittsburgh. PA 15206 | U.S.

Hardware	Qty. Needed	Source	ce	Pkg. Size	Pkgs. To Order	ost Per Unit	Total Cost	Parts for which Hardware is used	Finish
0.75in Steel Rod	1	McMaster Carr	<u>1346K350</u>	1	1	\$ 57.99	\$ 57.99		Steel
0.75in Sleeve Bearing	14	McMaster Carr	2938T540	1	14	\$ 1.39	\$ 19.46	Steel Rod	Bronze
0.75in Clamping Shaft Collar	2	McMaster Carr	6435K160	1	2	\$ 2.89	\$ 5.78	Steel Rod	Black-Oxide
3in Caster	4	McMaster Carr	2406T340	1	4	\$ 7.30	\$ 29.20		
5.50in Oval Grip Handle	2	McMaster Carr	1646A220	1	2	\$ 2.74	\$ 5.48		Dull Black
10-24 x 2.25in Hex Head Screw	8	McMaster Carr	91253A256	25	1	\$ 9.82	\$ 9.82	Handles	Black-Oxide
10-24 Locknut	8	McMaster Carr	90631A011	100	1	\$ 3.31	\$ 3.31	Handles	Zinc-Plated Steel
#10 Washers	8	McMaster Carr	98023A114	100	1	\$ 7.86	\$ 7.86	Handles	Zinc Yellow- Chromate Plated Steel
1/4-20 x 3.75in Hex Head Screw	16	McMaster Carr	91251A087	10	2	\$ 11.42	\$ 22.84	Casters	Black-Oxide
1/4-20 Locknut	16	McMaster Carr	90640A129	100	1	\$ 4.27	\$ 4.27	Casters	Zinc-Plated Steel
1/4in Washer	32	McMaster Carr	98023A029	100	1	\$ 7.70	\$ 7.70	Casters	Zinc Yellow- Chromate Plated Steel
400lb Tie-Down Ring	6	McMaster Carr	3076T340	1	6	\$ 1.50	\$ 9.00		Zinc-Plated Steel
#14 x 2in Wood Screw	12	McMaster Carr	90031A319	100	1	\$ 11.56	\$ 11.56	Tie-Down Rings	Zinc-Plated Steel









### International Society of Wheelchair Professionals I wheelchairnet.ORG 6425 Penn Avenue. Suite 401. Pittsburgh. PA 15206 | U.S.

Strap Hinge	3	McMaster Carr	<u>1530A530</u>	1	3	\$ 5.41	\$ 16.23		Dull Steel
#12 x 0.75in Wood Screw	24	McMaster Carr	90031A294	100	1	\$ 8.04	\$ 8.04	Hinges	Zinc-Plated Steel
10-24 x 5" Screw	2	McMaster Carr	91253A468	10	1	\$ 11.32	\$ 11.32	Stop bar	
#8 x 3in Construction Screw 10lb Box	200	Home Depot	3GCS10BK	722	1	\$ 32.47	\$ 32.47	Frame	Gold

McMaster Carr Total \$ 229.86 Home Depot Total \$ 32.47 Total Hardware Cost \$ 262.33

Lumber	Qty. Needed	Sourc	e	Pkg. Size	Pkgs. To Order	Cost Per Unit	Total Cost	Parts for which Hardware is used	Finish
2in x 6in x 12ft Lumber	10	Home Depot	<u>161756</u>	1	10	\$ 8.10	\$ 81.00	Frames	
4in x 6in x 8ft Lumber	1	Home Depot	<u>259270</u>	1	1	\$ 14.76	\$ 14.76	Frames	
0.75in x 4ft x 4ft Plywood	3	Home Depot	<u>454559</u>	1	3	\$ 45.98	\$ 137.94	Ramp and Top Sheet	Sande

Home Depot Total \$ 233.70 Total Lumber Cost \$ 233.70

Total Cost \$ 496.03










# **Drawings**







ITEM NO.	PART NUMBER	QTY.
1	Part A-Outside Left	1
2	Part A-Outside Right	1
3	Part A-Inner	2
4	Part B-Front	1
5	Part B-Back	1
6	Part C	10
7	Part C-With rib	3
8	Part D	2
9	Part E	1
10	Part F	3
11	Part G	3
12	Part G-Without wheel mount	1





(International):

(https://creativecommons.org/licenses/by-sa/4.0/)


FINISH: REV: DATE: 12/13/2019 1:32 Α WEIGHT (LBS): SHEET 1 OF 1 DO NOT SCALE DRAWING

В

Α

DESCRIPTION	PART NUMBER	QTY.
3/4" Steel rod	1346K350	1
3/4" Sleeeve bearing	2938T540	14
3/4" Shaft collar	6435K160	2
3" caster	2406T340	4
1/4"-20x 3-3/4" Screw	91251A087	16
1/4" Washer	98023A029	32
1/4"-20 Locknut	90640A129	16
Pull handle	1646A220	2
10-24 x 2-1/4" Screw	91253A256	8
10-24 Locknut	90631A011	8
10# Washers	98023A114	8
	3/4" Steel rod 3/4" Sleeeve bearing 3/4" Shaft collar 3" caster 1/4"-20x 3-3/4" Screw 1/4" Washer 1/4"-20 Locknut Pull handle 10-24 x 2-1/4" Screw 10-24 Locknut	3/4" Steel rod 1346K350 3/4" Sleeeve bearing 2938T540 3/4" Shaft collar 6435K160 3" caster 2406T340 1/4"-20x 3-3/4" 91251A087 Screw 98023A029 1/4"-20 Locknut 90640A129 Pull handle 1646A220 10-24 x 2-1/4" 91253A256 10-24 Locknut 90631A011

All screws used to put lumber together are #8 x3in Screws unless otherwise specified in the Bill of materials.







Copyright 2019, University of Pittsburgh. Made available under Creative Commons Attribution-ShareAlike 4.0 License (International):

(https://creativecommons.org/licenses/by-sa/4.0/)

TITLE:

#### Frame Hardware-Bill of Materials

UNLESS OTHERWI	SE SPECIFIED: DIMENSIONS ARE IN INCHES [mm]
MATERIAL:	FINISH:

TOLERANCES:	+	
SCALE:	SIZE:	DATE:

1:32	Α	12/13/2019
DO NOT SCALE DRAWING		WEIGHT (LBS):

.

REV:

SHEET 1 OF 1


В


2

В

Α

ITEM NO.	PART NUMBER	QTY.
1	Part J	1
2	Part K	1






Attribution-ShareAlike 4.0 License (International): (https://creativecommons.org/licenses/by-sa/4.0/) Α

В

TOLERANCES:	-	-	
SCALE:	SIZE:	DATE:	REV:
1:16	Α	12/13/2019	
DO NOT SCALE DRA	WING	WEIGHT (LBS):	SHEET 1 OF 1

FINISH:

ITEM NO.	PART NUMBER	QTY.
1	Part L	2
2	Part M	2
3	10-24 x 5" Screw. Part Number:91253A468	2



ISWP SM International Society of Wheelchair Professionals





Copyright 2019, University of Pittsburgh. Made available under Creative Commons Attribution-ShareAlike 4.0 License (International): (https://creativecommons.org/licenses/by-sa/4.0/) TITLE:

#### Stop bar- Bill of Materials

В

Α

REV:

SHEET 1 OF 1

UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES [mm]

MATERIAL: FINISH:


TOLERANCES: +


DO NOT SCALE DRAWING

SCALE:	SIZE:	DATE:	
1:12	Α		12/13/

12/13/2019
WEIGHT (LBS):

ITEM NO.	PART NUMBER	QTY.
1	Part H	1
2	Part I	1
3	3076T340_TIE-DOWN RING	6





В

Α

